

TimsPy: access timsTOF Pro data easily from Python



Mateusz K. Łacki¹, Sven Brehmer², Ute Distler¹, Stefan Tenzer¹


¹University Medical Center of the Johannes Gutenberg University Mainz, ²Bruker Daltonic GmbH

- timsTOF Pro is a novel instrument by Bruker
- uses trapped ion mobility spectrometry (TIMS) coupled with liquid chromatography (LC) and mass spectrometry (MS)



- collected data is stored in vendor's format
 - ★ accessible with freely available software development kit from many programming languages

- data dimensions include:



- but are recorded and stored as these numbers.

- TimsPy** gives you simple access to this data, the way data-scientists love it:

```
from timsipy import TimsDIA
D = TimsDIA('path/to/data_folder.d')
D[ frames, scans ]
```

```
Above, frames In: D[1:100,0:918]
select some 0 1 33 312260 9
frames and scans 1 1 34 220720 9
some scans, like 2 1 34 261438 9
on the right: 8394 99 914 313598 9
8395 99 917 354548 9
[1626073 rows x 4 columns]
```

Of course, we are talking of gigabytes of data. To minimize RAM usage, **TimsPy** offers iterators: `it = D.iter[1:10001, 0:918]`

```
In: next(it)
   frame scan tof_idx i
0 1 33 312260 9
1 1 34 220720 9
2 1 34 261438 9
...
1599 1 915 233374 9
1600 1 916 335348 77
[1601 rows x 4 columns]

In: next(it)
   frame scan tof_idx i
0 2 33 97298 9
1 2 33 310524 9
2 2 34 127985 9
...
6596 2 913 56442 9
6597 2 915 172202 9
[6598 rows x 4 columns]
```

and so on, until 10000th frame is reached

frames and scans can also be very general expressions, covering broad use-cases

```
D[1:5, [33, 50]]
D[(i**2 for i in range(10)), [33, 50]]
D[[1,2,10], 10:100]
D[[1,2,10], [10, 50]]
D['rt > 10 and rt < 50', 1:599]
```

And if you want to directly use physical quantities in the query?

We have you covered. Use:

```
D.phys[ RT, IM ] or
D.physIter[ RT, IM ]
```

```
In: D.phys[0.2:10]
   rt im mz
0 0.297934 1.601142 302.347671
1 0.297934 1.601142 1165.327281
2 0.297934 1.600000 391.984100
...
8936 9.995394 0.603017 266.728249
8937 9.995394 0.603017 1106.810657
[1554410 rows x 3 columns]
```

We still experiment with usage of VAEX and HDF5 to optimize your daily experience with timsTOF data on your laptop. You don't need a server or a supercomputer to study you data. Follow our [github](#) page!

- Give it a go now! It's as simple as:** `pip install timsipy`